



Manufacturers of pH & D.O.  
Sensors  
for Science and Industry


19 Thomas  
Irvine, CA 92618 • USA

Toll-Free:  
1-800-288-2833  
Tel: (949) 829-5555  
FAX: (949) 829-5560

## 12 mm OxyProbe Membrane Instructions


### STEP 1

While holding the membrane cartridge upright (with opening up) in one hand, carefully pour the electrolyte solution from the bottle into the cartridge until level is approximately 1/8" (3 mm) from the top edge of the cartridge. Remove any entrapped air bubbles by gently tapping on the side of the membrane cartridge sleeve. Try not to spill any of the electrolyte solution and be certain to wipe off any excess electrolyte from the outside surface of the cartridge sleeve with a clean tissue.



### STEP 2

In one hand, grip the membrane cartridge, which has been properly filled with electrolyte solution, by the thumb and fingers while holding the sensor assembly with the other hand as shown.



### STEP 3

Gently slide the membrane cartridge over the cathode / anode assembly until the cartridge threads seat against the mating threads of the sensor assembly. Thread the membrane cartridge onto the sensor assembly until there is a flush fitting between the cartridge sleeve and the tube. With a proper fit, the sealing o-ring should no longer be visible.

**Note:** Carefully inspect o-rings for any physical damage or excessive wear. Replace as needed or if there is any doubt as to their condition.